



# Transmitter implantation

S: OSI Part Number A3048S2-AA-C45-D

B: OSI Part Number A3049B2-AA-B45-B



- > Expose the skull
- > Clear of connective tissue
- > Clean (saline) and dry skull
- > Level skull

- > Key the skull with a needle.
- > Glue the skin to frontal bone
- > Fuse the frontal bones together\*
- > Drill holes (AP -3, ML +/-3)
- > Insert transmitter
- > Feed wires through neck
- > Straighten the coiled wire to form a hook long enough to touch cortex

- > Insert cannula into the drilled holes whilst pinning the wires into the hole.
- > Affix the cannulas to the skull with Relyx 3M (blue light)
- > Affix the wires to the skull with Medbond

- > Form the headplate using simplex cement
- > Screw on cannula caps

# Know your glue properties

GLUture – KA surgery wound closing



- >Flexible hold that dissolves in 5 days
- >Great for wound closing on scalp
- >NOT suitable for headplates or anything permanent

Medbond – sealing skin and silicon to skull and bone to bone



- >Not flexible, slow to dissolve
- >Great for sealing skin to bone
- >Excellent for bonding silicone of transmitter wires to the skull
- >Should be keyed after applying
- >Is suitable for base of headplates
- >If used in excess on skin can cause inflammation

3M Relyx Unicem – Adhering plastic cannula to bone



- >Rapid curing under blue light
- >Ideal for adhering cannula in place over bone.
- >Probably not as strong a bond as simplex
- >Does not bind well to silicon of transmitter wires.
- >Used in excess the transmitter wires will slip out the headpiece!

# Transmitter implantation

- > The A3049 transmitter requires a large pocket on the flank of the animal. Animals can access this site, overgroom and the skin in this area has a habit of sticking/putting pressure to the transmitter causing loss of blood supply and necrosis.
- > The A3048 transmitter is best situated along the spine.
- > Check cannula sizes

